Stiftung Tierärztliche Hochschule Hannover University of Veterinary Medicine Hannover, Foundation



EAAP Annual Meeting 2018



Sound Hooves:

Detection of lameness in dairy cows by acoustic analysis

N. Volkmann¹, B. Kulig², N. Kemper¹

¹Institute for Animal Hygiene, Animal Welfare and Farm Animal Behaviour ²Section of Agricultural and Biosystems Engineering; University of Kassel

Dubrovnik, 29.08.2018



Introduction



- Detection of lameness = important indicator for animal welfare!
- Currently:

No practical systems
for an early and
automated diagnosis
of lameness!





Introduction




Claw lesion

- → Pain
- → Altered weight distribution
- → Shorter and careful steps

(NEVEUX et al., 2006)

→ Systems working with weighing platforms or pressure sensitive walkways

(PASTELL et al., 2010; MAERTENS et al., 2011)






Aim



- Each step on solid surface produces a sound signal
- Signal differs between healthy/unhealthy claws



1. Duration:

- Lame cows do not run!
- 2. Maximum Volumen:
- Lame cows do not do such powerful steps

University of Veterinary Medicine Hannover, Foundation





- 77 lactating cows

Group 1: slatted floor, milking robot

Group 2: deep straw, milking parlor



- Alley with parts of slatted floor
- Cameras for sideview and on the limbs
- Animal weigh scale








Piezoelectric Sensors:

 Measurement of changes in pressure, strain or vibration and converting them to an electrical

charge

In this case:Recordingsound pressure (dB) !








Locomotion scoring

- LS0 = non-lame
- LS1 = lame

Claw lesions at hoof trimming

- Groups of diseases
 - 0 = none
 - 1 = non-infectious
 - 2 = infectious
 - 3 = both (infectious and non-infectious)







Parameters:

- 1. Duration of walking on the section of measuring (WS)
- 2. Standard deviation of volume in the recorded signal (StDevLoud)

Statistics:

SAS 9.4 (SAS INSTITUTE, 1999)
Comparison of means (PROC TTEST; PROC GLM)























Variable	mean	min	max	STD
WS (s)	3.71	2.12	5.82	0.78
StDevLoud (dB)	0.019	0.009	0.049	0.007


Analysed factors (P-values):

	Locomotion Score	Group of diseases	Type of pen	Weight
WS	< 0.0001	0.0812	0.5863	0.1829
StDevLoud	0.0010	0.0430	0.3120	0.5711






- Mean WS was faster in cows with LSO (nonlame) compared with score LS1 (lame)
- Conformation of LS:
 Cows with a smoother and quicker gait pattern got the lower score (FLOWER et al., 2006)







 Significant difference of volume (STDevLoud) between healthy animals (Group 0) and those with non-infectious diseases (Group 1)



 Cows with non-infectious diseases have a greater sensitivity to pain (PASSOS et al., 2017)



Conclusions



- Walking speed allows to draw conclusions to a smooth gait pattern. But: caution is required when using speed for lameness detection!
- If standard deviation of volume in the recorded signal (StDevLoud) can be interpreted as a variable for weight load of footsteps...
- ...first results clearly show the potential of walking sound analysis for lameness detection.



Outlook



- Recordings in the regular environment of the cows
- Testing of other factors (e.g. pollution of the claws)
- Clearer description of the sound signal
 - Number of steps
 - Stride duration and length
- Allocation of the sound to the individual limb
- Development of an individual (healthy) pattern for each cow





Thank you for your attention!

With support from





by decision of the German Bundestag

The project is supported by funds of the Federal Ministry of Food and Agriculture (BMEL) based on a decision of the Parliament of the Federal Republic of Germany via the federal Office for Agriculture and Food (BLE) under the innovation support programme.













